



Quality from Germany



User-friendly construction and high flexibility combined with maximum comfort make this blood donation/treatment chair optimal for blood plasma and whole blood donations, as well as for the therapy sector – not only for operators, but also for donors and patients.

All Likamed products are made in Germany and received certification marks like “TÜV” for top quality. Certification according to ISO 9001 and EN ISO 13485 is the appreciation for innovative and useful developments and for an outstanding product quality.



Designed for you, with you. The blood donor/treatment chair has been specially developed in line with the wishes and ideas of operators and donors.

- Infinitely adjustable chair surface by means of a 24V electric drive.
- Assumes the shock position quickly (5 sec.).
- Optimum positioning of the arm due to wedge-shaped arm rest.
- Arm rests can be rotated, and are adjustable in height.
- Upholstery and frame are easy to clean and disinfect.
- Changeable upholstery.
- Paper-roll holder
- Mobility thanks to large centrally locking castors.
- The upholstery is resistant to disinfectant and flame retardant. Ultra-high frequency welded and damp-proof upholstery quilting.

Options: (please ask for our options brochure)


- Fold down arm rest
- rechargeable battery pack
- Terry protective covers for cushions and armrests.
- Transparent plastic covers for arm-rest and foot rest.

Standard colours:



- Upholstery medium grey with matching grey arm and foot rests. Upholstery pastel blue with matching dark blue arm and foot rests. Frame and panelling in grey colour.


Other upholstery colours:



- red, yellow, mint, combined with grey or blue arm and foot rests. Frame and panelling in grey colour.



Individual adjustability of the arm rests and central locking (*=option).



Technical details and Guarantee

Armrests width	15 cm
Armrest length	53 cm
Laying surface width	55 cm
Total length	177 cm
Laying surface	188 cm
Total width with armrest	82 cm
Weight	48 kg
Access height	62 cm
Padding and neck support upholstery:	imitation leather flame retardant
Frame colour	lightgrey, RAL 7047
4 castors with locking brake, 1 x anti-static	central lock, Ø 10 cm
Transformer protection class	2, Typ BF
Transformer rated output	primary 230V/50 Hz, 110V/50 Hz

Transformer rated output	Secondary	25 V
Transformer rated output	Fuses	1 AT
ON duration		approx 10 %
EMV (electrical safety)	DIN EN 60601/1, 60601/1-2	
Protection class	II	
protection category	IPX4 (splash-proof)	
Drives	1 x 24 volt with limit circuit breaker	
Maximum load	150 kg	
Relative air humidity	to 30 - 75 %	
Ambient temperature	3 - 40°C	
Mains cable rolled	2,5 m	
Guarantee	2 Years	