



THE CHAIR FOR GREATEST DEMANDS



Seat part separately adjustable!

**200kg  
MAXIMUM LOAD**



Quality from Germany



**NOVO XS 3/4** – the next logical step in the expansion of our NOVO X model series. The NOVO X series not only offers exceptional quality, but also combines innovative technology with elegant design – this makes the entire range so very special. Thanks to the overall design, which has been intelligently developed to meet practical requirements; this product guarantees outstanding comfort for the patient as well as for the doctor and nurse.

**NOVO XS 4** – the top model as fully equipped treatment chair, leaving nothing to be desired. Four electric drives for separate, step-free adjustment of the back part, foot part, and at this model, even of the seat part, as well as a separate, step-free height adjustment of the complete laying surface. All movements are individually controllable by a hand switch. Easy to operate manual control switch with operate locking device.

**NOVO XS 3** – with all features of the XS 4 model, but without height adjustment. Both models are fitted with a flame retardant, resistant to disinfectant and toxic tested imitation leather upholstery.

All Likamed products are made in Germany and received certification marks like "TÜV" for top quality. Certification according to ISO 9001 and EN ISO 13485 is the appreciation we deserve for innovative and useful developments and for an outstanding product quality.

*For further information about the other features please take a look at our NOVO X brochure.*





## NOVO XS 4 – customised model for the recovery ward

This customised model was specially designed for the recovery ward. Alongside the tried and tested features seen in the XS4, this model also offers additional well-designed and extensively tested extras. And of course, the chair has been constructed to allow patient transfer.

### Additional features of this customised model for the recovery ward

(compared to NOVO XS 3/4)

- Gripping handle on backrest and NOVO XS 3/4 foam strip on the foot rest (with edge protection to prevent injuries) for operating the patient transfer switch.
- Battery pack
- Bracket for infusion apparatus stand.
- Infusion apparatus stand.
- The foot rest can be extended an additional 10 cm compared to the NOVO X model series.
- In the bed position the armrests are parallel to and level with the cushion.
- Central locking
- Ripping handle on the backrest for transport.

Other upholstery colours can be found in our options brochure.

Technical details and Guarantee NOVO XS 3/4 and NOVO XS 4 – customised model\*



Hessingpark Clinic Augsburg



Model	NOVO XS 4 customized*	NOVO XS 3/4
Armrests width		15 cm
Armrest length		55 cm
Laying surface width		60 cm
Total length	223 cm*	213 cm
Total length with extended footrest	230 cm*	220 cm
Width between armrests		60 cm
Sitting- laying height NOVO XS 4		54 - 80 cm
Sitting- laying height NOVO XS 3		58 cm
Weight		80 kg
Padding and neck support upholstery	Skai, flame retardant	
Colour: upholstery / neckrest	middle grey or pastel blue	
Colour: armrest / footrest	dark grey or dark blue	
Frame colour	light grey, RAL 7047	
4 castors with locking break		Ø 10 cm
Transformer	protection class II, Type BF	
Transformer rated output primary	230V/50 Hz, 110V/50 Hz	
	Secondary	25 V

Transformer rated output Fuses	1 A
ON duration	approx 10 %
EMV (electrical safety)	DIN EN 60601/1, 60601/1-2
NOVO XS 4	4 x 24 volt drive with limit circuit breaker
NOVO XS 3	3 x 24 volt drive with limit circuit breaker
<b>Maximum load</b>	<b>200 kg</b>
Relative air humidity	to 30 - 75 %
Ambient temperature	3 - 40°C
Mains cable rolled	2 m
Mains cable stretched	3,5 m
Back section tilting angle	up to 75°
Foot section tilting angle	up to 35°
protection category	IPX4 (splash-proof)
Protection class II	Guarantee
2 Years	
equipotential bonding	
Protection class type CF	optional