



NOVO X 550 – THE LIKAMED DIALYSIS-BED-CHAIR



NEW! Electrical adjustable foot part!



Quality from Germany



The all rounder: This, special for the dialysis developed bed-chair is an alternative to a bed or a recliner. A hand switch with a shock button and blind-friendly symbols allows the patients, thanks to five electric power systems, multifaceted, individual adjustability. For the safety all functions can be locked with a magnet key.

All Likamed products are made in Germany and received certification marks like "TÜV" for top quality. Certification according to ISO 9001 and EN ISO 13485 is the appreciation for innovative and useful developments and for an outstanding product quality.



OUR CUSTOM MADE DIALYSIS-BED-CHAIR
FOR MAXIMUM PATIENT COMFORT AND ULTIMATE SAFETY

Technology

- The new adjustable foot part allows to move into a sitting position.
- Step-free electromotive (5 electrical drives) adjustment of high stroke, back-, seat-, foot part and footrest.
- Rotating- and height adjustable armrests, alternatively left and right mountable.
- Four easy going, high-amperage double-wheels, central lockable.
- Robust steel frame with all the benefits of the LIKAMED range of loungers.
- Easy removable headpart
- Storage battery for network-independent operation.
- Hand switch with shock button, blind-friendly symbols and magnet key to lock the complete functions.



Comfort

- A low access height of 53 cm relieves the get on and about for shorter patients as well as the conversion of wheel-chair patients.
- An 40 cm large high stroke, guarantees an optimal working height for the operator.
- Extreme comfortable foam mattress. By request is a latex mattress as optional equipment available.
- Beneficial, space saving measurements, (fits through normal doorways).



Hygiene

- Subframe easily accessible.
- All parts easy to clean and disinfect.
- Strong long-lasting polyurethane armrest with hygienic surface.
- PU cover easy to clean.

Options

A range of carefully designed extras is available for the dialysis bed chair:

- Armrests with adjustable tilt
- Battery pack (for off-line use)
- Latex cover (mattress), 12 cm thick, CFC-free foamed
Cover: PU with zips across and down
- Reading lamp
- Terry stretch cover
- 2nd armrest
- Side bars
- Sit up help

Technical data

Base width	76 cm
Base length	213 cm
Width of bed (excl. leg and armrests)	84 cm
Overall width	94 cm
base height, incl. mattress (11 cm)	53-93 cm
Colour of frame	lightgrey RAL 7047
4 castors, central locking	Ø = 10 cm
5 electric motors with limit switch	24 V, IP 54
Transformer nominal power	100 VA, IP 54
Primary	230 V/ 50/ 60 Hz
optional	110 V/ 50/ 60 Hz
Fuse	1 Amp, slow-blow
Secondary	25 V
Handscharter	IP 54
Maximum load	175 kg
Weight	98 kg
Base adjustment:	
Back rest tilt	to 75°
Leg rest tilt	to 35°
Room temperature	3 - 40°C
Relative humidity	30 - 75 %

Art-Nr.: 13-00-01 / Index E 10/07