



## SMART - THE COMPACT ALL-ROUND ARM CHAIR FOR COUNTLESS MEDICAL USES



Quality from Germany



Choose the inexpensive solution! The compact construction and uncomplicated operation of the SMART medical arm chair will convince you that **SMART** is the best choice for a wide variety of medical use, an intelligent combination of modern technical performance unbeatable value for money.

All Likamed products are made in Germany and received certification marks like "TÜV" for top quality. Certification according to ISO 9001 and EN ISO 13485 is the appreciation we deserve for innovative and useful developments and for an outstanding product quality.



## SMART - THE ARM CHAIR THAT ENABLES YOU TO WORK EFFICIENTLY AND ECONOMICALLY


- Seat and base progressively and electrically adjustable.
- Extensive adjustable from normal sitting position to shock-position.
- Seat/base surface designed and shaped ergonomically.
- Tilting polyurethane foam armrests, height adjustable.
- Easy to clean imitation leather, hardly inflammable.
- Hygienical surfaces on all parts, easy to clean and disinfect.
- Frame in modern, friendly light grey RAL 7047.

### Standard colours:




- Upholstery middle grey with matching dark grey arm- and footrests.
- Upholstery pastel blue with matching dark blue arm- and footrests.
- Frame and panelling in grey colour.

### Other upholstery colours:



- red, yellow, mint, combined with grey or blue arm and foot rests.



### Options:

Please ask for our special brochure.

### Technical details:

Armrest width	12 cm
Armrest length	52 cm
Laying surface width	55 cm
Total length	188 cm
Total width	77 cm
Weight	48 kg
Padding and neck support upholstery:	imitation leather flame retardant
Colour: Footpart and armrest	standard dark grey
Colour upholstery:	middle grey
Frame colour	light grey RAL 7047
four castors, all lockable	Ø 7,5 cm
Transformer rated output	Primary 230 V / 50/60 HZ

Transformer rated output	Primary	110 V / 50/60HZ
Transformer rated output	Secondary	25 V
	Fuse	1 A
ON duration	approx. 10 %	
Protection class	II	
EMC (electrical safety)	DIN EN 60601/1, 60601/1-2	
Drive	1 x 24 V with limit switch	
Maximum load	150 kg	
Relative air humidity	t o 30 - 75 %	
Ambient temperature	3 - 40°C	
Main cable stretched	3 m	
Guarantee	2 Years	