

LIKAMED vaccination chair

The compact chair for
your vaccinate patients

Make sure
to get the
price bonus
of up to
EUR 200



Optimally equipped

The compact therapy chair in its specially developed configuration for vaccinations offers price- and performance-oriented comfort in TÜV-tested quality.

Stay well equipped for the efficient treatment and safe positioning of your patients due to the infinitely variable motorised adjustment up to the shock position, combined with mobility and stability, even in patients with high body weight. With its integrated paper roll holder, it provides you with a practical hygienic solution for a quick patient turnaround.



The vaccination chair is ideally suited for:

- Vaccination centres
- General medical practices
- Pharmacies
- Nursing homes
- Blood and plasma donation
- and many other treatment areas

Accessories for the complete
vaccination station:

- Work and patient table
- Work stool



likamed.de/vaccination-chair

find out more



LIKAMED®

Main features of the LiKAMED vaccination chair



Flexible: sitting, comfort and shock position

The infinitely variable motorised adjustment of the reclining surface lets you bring the patient directly into the shock position in an emergency.



Convenient: hand switch

For infinitely variable adjustment of the reclining surface. With manual key lock function and fastening clip on the back.




Mobile: Casters Ø 7.5 cm with individual lock

The low access height simplifies entry and exit for the patient and assistance by the nursing staff. The individual lock on all four casters keeps it stationary.



Fast turnaround of patients: paper roll holder included

Hygiene and cleanliness for every patient by changing the paper cover. The permanently installed roll holder keeps the medical paper roll ready for use and positions the paper cover.



Additional protection: antimicrobial upholstery included

The special antimicrobial imitation leather can reduce the viability or infectiousness of microorganisms, thus offering greater safety for patient and health care personnel. The upholstery is available in seven colour options.



Relaxed patients: neck roll included

The neck roll in upholstery colour provides relaxation in the neck area.

more information: likamed.de/vaccination-chair

Complete your vaccination workstation:



Work stool

- Full foam seat part or foam seat part covered with imitation leather
- Seat part with easy ring control
- Infinitely adjustable via gas spring
- Soft casters for hard floors
- 6 colours for foam seat part



Work and patient table

- single-handed height adjustment
- large removable tray with recessed grips
- with integrated disinfectant bottle holder in the tray
- removable plastic waste container including lid
- hanger for bags
- **integrated bottle holder included**



 CONTACT PERSONS

 CONTACT FORM

 E-MAIL

LIKAMED[®]

LIKAMED GmbH · Raußmühlstraße 7 · 75031 Eppingen

Phone: +49 7262 9189-0 · Fax: +49 7262 9189-900 · E-mail: info@likamed.de · www.likamed.de/en