

SENSA® i Scale

Optimally equipped
for therapy



New
improved
model!

LIKAMED®

SENSA® i Scale

Therapy chair with integrated scale



The calibrated or certified calibrated scale is suitable for a safe working load of up to 200 kg and the four load cells enable weighing of the patient in the preferred position, considerably easing the work of healthcare personnel.

Integrated scale

- SENSA i Scale: calibrated scale
- SENSA i Scale+: certified calibrated scale (required in the EU)
- 4 double-beam load cells with strain gauge application (load cells C3)
- Scale range 2–200 kg (resolution 100 g)
- Integrated spirit level in the base frame to ensure correct positioning

Display unit

- With special dialysis function button
- Illuminated 7-segment display
- Digit height 20 mm
- With RS232 interface

Classification as per 2017/745 (MDR): Im (with measuring function)



tuv-sud.com/ps-cert

CE 0123



Made in Germany



Ease of operation

Hand switch with comfort and special buttons and magnetic key lock function.



Flexible features

Configure your SENSA® i Scale

- Central locking
- Infusion rod holder
- Electrically adjustable footrest
- and much more



Lots of new options!

scale.configurator.likamed.de

© LIKAMED GmbH, 61-PR-1527-EN 2021.11

The information is up to date at the time of printing. LIKAMED GmbH accepts no liability for differences from the illustrations regarding colour and format. No responsibility is accepted for errors and misprints. LIKAMED GmbH reserves the right to make changes to publications and products at any time. Where the masculine form or the generic masculine is used in this brochure for the sake of legibility, this specifically includes the feminine and other gender identities.

LIKAMED®

LIKAMED GmbH
Raußmühlstraße 7 · 75031 Eppingen · Germany
Phone: +49 7262 9189-0 · E-mail: info@likamed.de
[likamed.de](https://www.likamed.de)

