


### Unique seating comfort

The new pressure relieving gel seat cushion by LiKAMED: the gel seat cushion enables with its high restoring force and the good thermal conductivity a cooling character and ensures a unique seating comfort. Due to the 3D deformability of the polyurethane-based gel, the pillow provides a pressure relief over the whole surface and ensures your patients more comfort during long treatment sessions.



## Pressure relieving gel seat cushion



- polyurethane-based 3D gel
- biocompatible and odorless
- free of plasticizers or other volatile substances
- water-repellent PU-cover in grey, easily removable and wipeable with a damp cloth
- backpart with non-skid material
- dimensions: 41 x 41 x 5 cm
- fits all LiKAMED treatment chairs

**LiKAMED®**

LiKAMED GmbH · Raußmühlstraße 7 · 75031 Eppingen · Germany

Phone: +49 7262 9189-0 · Fax: +49 7262 9189-900 · E-Mail: [info@likamed.de](mailto:info@likamed.de) · [www.likamed.de](http://www.likamed.de)